
scikit-surgeryarucotracker Documentation

Stephen Thompson

Apr 27, 2023

Contents

1	scikit-surgeryarucotracker	3
1.1	Installing	4
1.2	Using	4
1.3	Developing	4
1.4	Licensing and copyright	5
1.5	Acknowledgements	5
2	ARuCo Tracking	7
3	Register detected image points (2D) to model points (3D)	9
4	Manage rigid bodies consisting of multiple tags	11
	Index	13

Source code is available on [GitHub](#).

CHAPTER 1

scikit-surgeryarucotracker

Author: Stephen Thompson

scikit-surgeryarucotracker provides a simple Python interface between OpenCV's ARuCo marker tracking libraries and other Python packages designed around scikit-surgerytrackers. It allows you to treat an object tracked using ARuCo markers in the same way as an object tracked using other tracking hardware (e.g. aruco - scikit-surgerynditracker).

scikit-surgeryarucotracker is part of the [SciKit-Surgery](#) software project, developed at the [Wellcome EPSRC Centre for Interventional and Surgical Sciences](#), part of [University College London \(UCL\)](#).

scikit-surgeryarucotracker is tested with Python 3.6 and may support other Python versions.

1.1 Installing

```
pip install scikit-surgeryarucotracker
```

1.2 Using

Configuration is done using Python libraries. Tracking data is returned in NumPy arrays.

```
from sksurgeryarucotracker.arucotracker import ArUcoTracker
config = {
    "video source" : 0
}
tracker = ArUcoTracker(config)

tracker.start_tracking()
print(tracker.get_frame())
tracker.stop_tracking()
tracker.close()
```

1.3 Developing

1.3.1 Cloning

You can clone the repository using the following command:

```
git clone https://github.com/SciKit-Surgery/scikit-surgeryarucotracker
```

1.3.2 Running the tests

You can run the unit tests by installing and running tox:

```
pip install tox
tox
```

1.3.3 Contributing

Please see the [contributing guidelines](#).

1.3.4 Useful links

- [Source code repository](#)
- [Documentation](#)

1.4 Licensing and copyright

Copyright 2019 University College London. scikit-surgeryarucotracker is released under the BSD-3 license. Please see the [license file](#) for details.

1.5 Acknowledgements

Supported by Wellcome and EPSRC.

A class for straightforward tracking with an ARuCo

class `sksurgeryarucotracker.arucotracker.ArUcoTracker` (*configuration*)

Bases: `sksurgerycore.baseclasses.tracker.SKSBASETracker`

Initialises and Configures the ArUco detector

Parameters configuration – A dictionary containing details of the tracker.

video source: defaults to 0

aruco dictionary: defaults to DICT_4X4_50

marker size: defaults to 50 mm

camera projection: defaults to None

camera distortion: defaults to None

smoothing buffer: specify a buffer over which to average the tracking, defaults to 1

rigid bodies: a list of rigid bodies to track, each body should have a 'name', a 'filename' where the tag geometry is defined, and an 'aruco dictionary' to use. Additionally we can include 'tag width' in mm when the tag has been scaled during printing or is displayed on a mobile phone screen or similar

Raises Exception – ImportError, ValueError

close ()

Closes the connection to the Tracker and deletes the tracker device.

Raises Exception – ValueError

get_frame (*frame=None*)

Gets a frame of tracking data from the Tracker device.

Parameters frame – an image to process, if None, we use the OpenCV video source.

Returns

`port_numbers`: If tools have been defined port numbers are the tool descriptions. Otherwise port numbers are the aruco tag ID prefixed with aruco

`time_stamps` : list of timestamps (cpu clock), one per tool

`frame_numbers` : list of framenumbers (tracker clock) one per tool

`tracking` : list of 4x4 tracking matrices, rotation and position,

`tracking_quality` : list the tracking quality, one per tool.

Raises Exception – ValueError

`get_tool_descriptions()`

Returns tool descriptions

`has_capture()`

:Returns true if the tracker has it's own opencv source otherwise false

`start_tracking()`

Tells the tracking device to start tracking. :raise Exception: ValueError

`stop_tracking()`

Tells the tracking devices to stop tracking. :raise Exception: ValueError

Register detected image points (2D) to model points (3D)

Classes and functions for 2D to 3D registration

`sksurgeryarucotracker.algorithms.registration_2d3d.estimate_poses_no_calibration` (*marker_corners*, *aruco_board*)

Returns tracking data for a camera with no calibration data. *x* and *y* are the screen pixel coordinates. *z* is based on the size of the tag in pixels, there is no rotation. No account is taken of the size of the model marker pattern, so it will be bit flakey.

`sksurgeryarucotracker.algorithms.registration_2d3d.estimate_poses_with_calibration` (*marker_corners*, *marker_ids*, *aruco_board*, *camera_projection_matrix*, *camera_distortion*)

Estimate the pose of a single tag or a multi-tag rigid body when the camera calibration is known. *marker_corners2d*: a list of 2d marker corners, 1 row per tag,

8 columns per tag

Parameters

- **model_points** – Matched list of of corresponding model points, 1 row per tag, 15 columns per tag: corner points and centre point
- **camera_projection_matrix** – a 3x3 camera projection matrix
- **camera_distortion** – camera distortion vector

return : a tracking rotation, translation and a quality

Manage rigid bodies consisting of multiple tags

Classes and functions for maintaining ArUco rigid bodies

class `skSURGERYarucotracker.algorithms.rigid_bodies.ArUcoRigidBody` (*rigid_body_name*)
 Bases: `object`

Class to handle the loading and registering of ArUco Rigid Bodies

add_single_tag (*tag_size, marker_id, dictionary*)

We can use this to track single ArUco tags rather than patterns as long as we know the tag size in mm

Param *tag_size* in mm

Param *marker_id* / _

get_dictionary_name ()

returns the name of the aruco dictionary in use

get_pose (*camera_projection_matrix, camera_distortion*)

Estimate the pose of the rigid body, with or without camera calibration

Param *camera_projection_matrix* 3x3 projection matrix. If None we estimate pose based on pattern size

Param 1x5 camera distortion vector

load_3d_points (*filename, dictionaryname*)

Loads the 3D point geometry from a file

Parameters filename – Path of file containing tag data

reset_2d_points ()

Clears 2D point lists.

scale_3d_tags (*measured_pattern_width*)

We can scale the tag, which is very useful if you've got the tag on your mobile phone.

Parameters measured_pattern_width – Width of the tag in mm

set_2d_points (*two_d_points, tag_ids*)

takes a list of 2 points, and if the id's match 3D points, assigns them to a list of 2d points

Parameters

- **two_d_points** – array of marker corners, 4 for each tag
- **tag_ids** – id for each tag

Returns tag ids for any assigned tags

class `sksurgeryarucotracker.algorithms.rigid_bodies.Board` (*markerpoints, dictionary, marker_ids*)

Bases: object

A local replacement for aruco.Board which was deprecated at 4.7

class `sksurgeryarucotracker.algorithms.rigid_bodies.TwoDTags`

Bases: object

Stores two linked arrays, one of tag IDs and the other 2D points

append_tag (*tag_id, points*)

Adds a tag to the two point list :param tag_id: The id of the tag :param points: 4 points defining the tag corners

`sksurgeryarucotracker.algorithms.rigid_bodies.configure_rigid_bodies` (*configuration*)
reads configuration and creates a list of rigid bodies together with a list of dictionaries used.

`sksurgeryarucotracker.algorithms.rigid_bodies.load_board_from_file` (*filename, dictionary=16*)

loads marker pattern from filename. :return: an aruco.board :raise ValueError: If the file does not have 16 or 13 columns

`sksurgeryarucotracker.algorithms.rigid_bodies.scale_tags` (*board, measured_pattern_width*)

We can scale the tag on a board, which is very useful if you've got the tag on your mobile phone.

Param the board to scale

Parameters **measured_pattern_width** – Width of the tag in mm

`sksurgeryarucotracker.algorithms.rigid_bodies.single_tag_board` (*tag_size, marker_id, dictionary=<cv2.aruco.Dictionary_0x7f8bac399930>*)

Create a board consisting of a single ArUco tag

Param tag size in mm

Param marker id

Param dictionary to use

A

`add_single_tag()` (*sksurgeryaruco-tracker.algorithms.rigid_bodies.ArUcoRigidBody method*), 11

`append_tag()` (*sksurgeryaruco-tracker.algorithms.rigid_bodies.TwoDTags method*), 12

`ArUcoRigidBody` (*class in sksurgeryaruco-tracker.algorithms.rigid_bodies*), 11

`ArUcoTracker` (*class in sksurgeryaruco-tracker.arucotracker*), 7

B

`Board` (*class in sksurgeryaruco-tracker.algorithms.rigid_bodies*), 12

C

`close()` (*sksurgeryaruco-tracker.arucotracker.ArUcoTracker method*), 7

`configure_rigid_bodies()` (*in module sksurgeryaruco-tracker.algorithms.rigid_bodies*), 12

E

`estimate_poses_no_calibration()` (*in module sksurgeryaruco-tracker.algorithms.registration_2d3d*), 9

`estimate_poses_with_calibration()` (*in module sksurgeryaruco-tracker.algorithms.registration_2d3d*), 9

G

`get_dictionary_name()` (*sksurgeryaruco-tracker.algorithms.rigid_bodies.ArUcoRigidBody method*), 11

`get_frame()` (*sksurgeryaruco-tracker.arucotracker.ArUcoTracker method*), 7

`get_pose()` (*sksurgeryaruco-tracker.algorithms.rigid_bodies.ArUcoRigidBody method*), 11

`get_tool_descriptions()` (*sksurgeryaruco-tracker.arucotracker.ArUcoTracker method*), 8

H

`has_capture()` (*sksurgeryaruco-tracker.arucotracker.ArUcoTracker method*), 8

L

`load_3d_points()` (*sksurgeryaruco-tracker.algorithms.rigid_bodies.ArUcoRigidBody method*), 11

`load_board_from_file()` (*in module sksurgeryaruco-tracker.algorithms.rigid_bodies*), 12

R

`reset_2d_points()` (*sksurgeryaruco-tracker.algorithms.rigid_bodies.ArUcoRigidBody method*), 11

S

`scale_3d_tags()` (*sksurgeryaruco-tracker.algorithms.rigid_bodies.ArUcoRigidBody method*), 11

`scale_tags()` (*in module sksurgeryaruco-tracker.algorithms.rigid_bodies*), 12

`set_2d_points()` (*sksurgeryaruco-tracker.algorithms.rigid_bodies.ArUcoRigidBody method*), 11

`single_tag_board()` (*in module sksurgeryaruco-tracker.algorithms.rigid_bodies*), 12

`sksurgeryaruco-tracker.algorithms.registration_2d3d` (*module*), 9

`sksurgeryaruco-tracker.algorithms.rigid_bodies` (*module*), 11

`sksurgeryarucotracker.arucotracker` (*module*), 7

`start_tracking()` (*sksurgeryarucotracker.arucotracker.ArUcoTracker* method), 8

`stop_tracking()` (*sksurgeryarucotracker.arucotracker.ArUcoTracker* method), 8

T

`TwoDTags` (*class* in *sksurgeryarucotracker.algorithms.rigid_bodies*), 12